

# ASD19-N

## Single Channel NIR Absorption Probe




- Inline real-time process monitoring
- Concentration measurement insensitive to color changes
- Superior sapphire window with no seals, gaps or crevices
- Extremely low maintenance
- CIP/SIP-compatible
- All wetted material and surface finishes are certifiable and traceable

The model ASD19-N is a high-precision, single channel absorption probe. The stainless steel probe is specially engineered for use in laboratory or pilot-scale fermenters or bioreactors. It precisely measures the cell growth of microbial or cell cultures as a function of NIR absorption. The extremely durable and ultra-sanitary ASD19-N probe is designed for autoclaving.

The ASD19-N probe uses the light in the Near Infrared range (NIR) from 840 to 910 nm. A precisely defined and constant LED light beam penetrates the process medium. The attenuation of the light intensity, caused by absorption and/or scattering by dissolved and undissolved substances, is detected by a hermetically sealed photodiode.

The ASD19-N probe is available with two different insertion depths and an optical path length (OPL) of 1 mm, 5 mm, 10 mm or 20 mm. Shorter optical path lengths are typically used for dense cell cultures such as bacterial and yeast cultures. Longer OPLs are used for lower absorbing mammalian cell cultures.

The seal-less, sapphire optical window design eliminates crevices and gaps to assure the highest level of sterility, cleanability and sensor integrity. The surface of the wetted parts is electro-polished and may optionally comply with BN2. The ASD19-N probe mounts conveniently through an industry-standard, 19 mm port (M26x1) in the fermenter or bioreactor head plate.



### Type ASD19-N

- 1 Optical path length (OPL)
- 2 Sapphire windows
- 3 Detector
- 4 Day-light filter
- 5 LED-Light source

# Technical Data

## ASD19-N Probe



**Material:**

Wetted parts: stainless steel 1.4435 (SS 316 L)  
 Surface: electro-polished Ra < 0.8 µm (standard)  
 Windows: sapphire (without gasket)  
 Housing: stainless steel 1.4571 (SS 316 Ti)

**Port gasket:**

O-ring Ø 15.60 x 1.78 mm

**Gasket material:**

application specific, selection by end user

**Permitted:**

EPDM (FDA), others on request

**Port connection:**

Fermenter cover plates  
 Diameter: Ø 19 mm (H7)  
 Shoulder: Ø 22 mm  
 Thread: M26x1

**Insertion depth maximal:**

- ASD19-N-EA: OPL + 210 mm from shoulder
- ASD19-N-EB: OPL + 310 mm from shoulder

**Optical path length (OPL):**

1, 5, 10 or 20 mm

**Pressure rating:**

PN10 (Test pressure PT 15 bar)

**Permitted pressure PS:**

10 mbar - 10 bar with TS 0 °C / +100 °C

**Permitted pressure at elevated temperature:**

|          |       |     |     |
|----------|-------|-----|-----|
| TS [°C]  | < 100 | 125 | 150 |
| PS [bar] | 10    | 8   | 6   |

**Permitted process temperature TS:**

- permanent: +5 °C to +65 °C, (+41 °F to +149 °F)
  - peak (60 min/day): +5 °C to +135 °C, (+41 °F to +275 °F)
  - peak (30 min/day): +5 °C to +145 °C, (+41 °F to +293 °F)
- autoclaving possible without cable, short term values are only valid if the LED is switched off

**Ambient temperature:**

- operation: 0 °C to +40 °C, (+32 °F to +104 °F)
- transportation: -20 °C to +70 °C, (-4 °F to +158 °F)

**Air purge:**

connections M5 available as standard

**Light source:**

Hybrid LED, hermetically sealed

**Wavelength range:**

840 - 910 nm

**Detector:**

silicon photodiode, hermetically sealed

**Measuring range:**

any measuring range between  
 0 - 4 CU

**Cable connection:**

probe cable ASD-CC, pluggable both sides  
 2, 3, 5 or 10 m, (7, 10, 16 or 33 ft.)

**Weight:**

- Probe: approx. 2.0 kg, depending on version

**Type of protection:**

IP65

**Certificates:**

ISO 9001:2000, PED, CE, HPO

**Use with Fermenter Control converter!**

## Options



**ASD19-N-EA-SR / ASD19-N-EB-SR**  
 wetted parts with surface electro-polished  
 Ra < 0.4 µm, dF < 1% (BN2)